ULTRASONIC TOOTHBRUSH

User manual
Declaration of warranty

Removing the rechargeable battery

WARNING


Only remove the rechargeable battery when you discard the appliance. Make sure the battery is completely empty when you remove it. To remove the rechargeable battery, you need a towel or cloth, a hammer and a flat-head (standard) screwdriver. Observe basic safety precautions when you follow the procedure outlined below. Be sure to protect your eyes, hands, fingers, and the surface on which you work. **Note:** We strongly advise you to have a professional remove the rechargeable battery.


b)


c)


d)


e), f)

g)

h)

i)

j)

- a) To deplete the rechargeable battery of any charge, remove the handle from the charger, turn on the electric toothbrush and let it run until it stops. Repeat this step until you can no longer turn on the electric toothbrush.
- b) Remove and discard the brush head. Cover the entire handle with a towel or cloth.
- c) Hold the top of the handle with one hand and strike the handle housing 0.5 inch above the bottom end. Strike firmly with a hammer on all 4 sides to eject the end cap. **Note:** You may have to hit on the end several times to break the internal snap connections.
- d) Remove the end cap from the toothbrush handle. If the end cap does not release easily from the housing, repeat step 3 until the end cap is released.
- e) Holding the handle upside down, press the shaft down on a hard surface. If the internal components do not easily release from the housing, repeat step 3 until the internal components are released.
- f) Remove the rubber battery cover.
- g) Wedge the screwdriver between the battery and the black frame at the bottom of the internal components. Then pry the screwdriver away from the battery to break the bottom of the black frame.
- h) Insert the screwdriver between the bottom of the battery and the black frame to break the metal tab connecting the battery to the green printed circuit board. This will release the bottom end of the battery from the frame.
- i) Grab the battery and pull it away from the internal components to break the second metal battery tab. **Caution:** Be aware of the sharp edges of the battery tabs so as to avoid injury to your fingers.
- j) Cover the battery contacts with tape to prevent any electrical short from residual battery charge. The rechargeable battery can now be recycled, and the rest of the product discarded appropriately.

Cleaning

TOOTHBRUSH HANDLE

- a) Remove the brush head and rinse the metal shaft area with warm water. Make sure you remove any residual toothpaste. **Caution:** Do not push on the rubber seal on the metal shaft with sharp objects, as this may cause damage.
- b) Wipe the entire surface of the handle with a damp cloth. **Note:** Do not use isopropyl rubbing alcohol, vinegar or bleach to clean handle as this may cause discoloration.

BRUSH HEAD

- a) Rinse the brush head and bristles after each use.
- b) Remove the brush head from the handle and rinse the brush head connection with warm water at least once a week.

Note: Do not clean the brush head in the dishwasher.

WIRELESS CHARGER BASE

- a) Unplug the charger before you clean it.
- b) Use a damp cloth to wipe the surface of the charger and travel case.

Storage

If you are not going to use your electric toothbrush for an extended period of time, unplug the charger from the wall socket, clean it and store it in a cool and dry place away from direct sunlight.

Safety warnings

**CAUTION!**

This symbol highlights dangerous situations which may lead to light injuries, or damage, malfunction, and/or destruction of the device.

Power supply

This product shall only be connected to an external power supply rated at 5V dc, and a minimum rated current of 400mA;

Any external power supply used with this product shall comply with relevant regulations and standards applicable in the country of intended use;

Connect the charging base to a laptop, a wall-mounted USB socket or a suitable wall adapter. The USB cord is attached to the charging base;

Important: The wall adapter must be a safe extra low voltage source, compatible with the electrical ratings marked on the device;

Don't leave your rechargeable battery discharged or unused for extended periods of time. To extend battery life, charge them at least every 6-9 months;

Charge the battery for above 12 hours before first use. Please understand that wireless charging is comparably slow for the battery charging;

Avoid storage at extreme temperatures (above 65°C/149°F or below -10°C/14°F). This can shorten battery life and may affect the product.

DANGER

Keep the charger away from water. Do not place or store it over or near water contained in a bathtub, washbasin, sink etc. Do not immerse the charger in water or any other liquid. After cleaning, make sure the charger is completely dry before you connect it to the mains.

WARNING

Check the USB cord before you connect the local mains;

Discontinue use if this product appears damaged in any way (brush head, handle, charger).

This product contains no serviceable parts.

Never use the charger if it has a damaged cord or plug. Keep away from heated surfaces. Do not use the charger outdoors;

Avoid using a brush head with crushed or bent bristles. Damaged bristles may break off while brushing. Replace the brush head with a new one every 3 months or sooner if signs of wear appear. Inappropriate cleaning may cause bristle damage;

This product is designed to clean your teeth, gums and tongue only. Use this product only for its intended use as described in this booklet. Discontinue use of this product and contact a physician/dentist if discomfort or pain is experienced;

This appliance is not intended for use by persons (including children) with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction concerning use of the appliance by a person responsible for their safety;

Children should be supervised to ensure that they do not play with the appliance;

This product is a personal care device and is not intended for use on multiple patients in a dental office or institution; Do not use attachments other than those recommended by the manufacturer;

If your toothpaste includes peroxide, baking soda or other bicarbonate (common in whitening toothpastes), thoroughly clean your brush head and handle with soap and water after each use. Otherwise plastic cracking may occur;

Do not clean the brush head or the charger in the dishwasher.

Medical warnings

Consult your dental professional before you use this product if you have had oral or gum surgery in the previous 2 months;
Contact your dental professional if excessive bleeding occurs after using this product or bleeding continues to occur after 1 week of use;
Address questions concerning usage with a pacemaker or other implanted device to your physician or the implant device manufacturer;
This Product has been tested and is compliant with safety standards for electromagnetic devices;
Consult your physician prior to using this product if you have medical concerns.

Customer Service Information


If you have any questions after reading this manual, feel free to contact our Customer Service Center or visit our website at www.prozis.com/support


SAVE THESE INSTRUCTIONS FOR FUTURE REFERENCE

EN


Quick Guide

Read the complete manual at www.prozis.com/user-manuals

EN

Statement of legal compliance

By means of this document, ONSALESIT declares that the Prozis Ultrasonic Toothbrush is fully compliant with all the standards specified in the applicable legislation as well as other relevant clauses.

The product's declaration of conformity is available on this address:

www.prozis.com/user-manuals.

For any additional explanation, ONSALESIT will always be available through our customer support communication channels, found at www.prozis.com/support.

To avoid any negative impact on the environment or on human health due to the dangerous substances present in electrical and electronic devices, the end users of such devices are expected to understand the meaning of the symbol consisting of an image of a crossed-out waste container. Do not dispose of waste electrical and electronic equipment along with unsorted household waste, instead dispose of it separately and properly.

Declaration of warranty

This device is covered by a limited warranty found in the corresponding instruction manual, that can be reviewed on the following URL: www.prozis.com/user-manuals.

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PROZIS / ONSALESIT, S.A.